



TRIM35 Antibody, FITC conjugated

Product Code	CSB-PA892166LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UPQ4
Immunogen	Recombinant Human Tripartite motif-containing protein 35 protein (1-206AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tripartite motif-containing protein 35 (Hemopoietic lineage switch protein 5), TRIM35, HLS5 KIAA1098
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TRIM35
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.