



PILRA Antibody, Biotin conjugated

Product Code	CSB-PA892331LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UKJ1
Immunogen	Recombinant Human Paired immunoglobulin-like type 2 receptor alpha protein (22-226AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Paired immunoglobulin-like type 2 receptor alpha (Cell surface receptor FDF03) (Inhibitory receptor PILR-alpha), PILRA
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	PILRA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.