




RARRES3 Antibody

Product Code	CSB-PA892340LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UL19
Immunogen	Recombinant Human Retinoic acid receptor responder protein 3 protein (1-133AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Retinoic acid receptor responder protein 3 (EC 3.1.1.-) (HRAS-like suppressor 4) (HRSL4) (RAR-responsive protein TIG3) (Retinoid-inducible gene 1 protein) (Tazarotene-induced gene 3 protein), RARRES3, RIG1 TIG3
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	RARRES3

Image



Immunohistochemistry of paraffin-embedded human lung cancer using CSB-PA892340LA01HU at dilution of 1:100