



# SLC7A7 Antibody, FITC conjugated

Product Code	CSB-PA892349LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UM01
Immunogen	Recombinant Human Y+L amino acid transporter 1 protein (325-382AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Y+L amino acid transporter 1 (Monocyte amino acid permease 2) (MOP-2) (Solute carrier family 7 member 7) (y(+)-L-type amino acid transporter 1) (Y+LAT1) (y+LAT-1), SLC7A7
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC7A7