



RPL26L1 Antibody, Biotin conjugated

Product Code	CSB-PA892361LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UNX3
Immunogen	Recombinant Human 60S ribosomal protein L26-like 1 protein (1-145AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	60S ribosomal protein L26-like 1 (Large ribosomal subunit protein uL24-like 1), RPL26L1, RPL26P1
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	RPL26L1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.