



SLC12A4 Antibody, Biotin conjugated

Product Code	CSB-PA892363LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UP95
Immunogen	Recombinant Human Solute carrier family 12 member 4 protein (601-900AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Solute carrier family 12 member 4 (Electroneutral potassium-chloride cotransporter 1) (Erythroid K-Cl cotransporter 1) (hKCC1), SLC12A4, KCC1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC12A4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.