



ACAD8 Antibody, HRP conjugated

Product Code	CSB-PA892457LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UKU7
Immunogen	Recombinant Human Isobutyryl-CoA dehydrogenase, mitochondrial protein (201-415AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Isobutyryl-CoA dehydrogenase, mitochondrial (EC 1.3.99.-) (Activator-recruited cofactor 42 kDa component) (ARC42) (Acyl-CoA dehydrogenase family member 8) (ACAD-8), ACAD8, ARC42 IBD
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ACAD8
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.