



SRPK3 Antibody, Biotin conjugated

Product Code	CSB-PA892486LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UPE1
Immunogen	Recombinant Human SRSF protein kinase 3 protein (1-300AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	SRSF protein kinase 3 (EC 2.7.11.1) (Muscle-specific serine kinase 1) (MSSK-1) (Serine/arginine-rich protein-specific kinase 3) (SR-protein-specific kinase 3) (Serine/threonine-protein kinase 23), SRPK3, MSSK1 STK23
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SRPK3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.