



POLDIP2 Antibody, HRP conjugated

Product Code	CSB-PA896496LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y2S7
Immunogen	Recombinant Human Polymerase delta-interacting protein 2 protein (1-200AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Polymerase delta-interacting protein 2 (38 kDa DNA polymerase delta interaction protein) (p38), POLDIP2, PDIP38 POLD4
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	POLDIP2