



POLDIP2 Antibody, Biotin conjugated

Product Code	CSB-PA896496LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y2S7
Immunogen	Recombinant Human Polymerase delta-interacting protein 2 protein (1-200AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Polymerase delta-interacting protein 2 (38 kDa DNA polymerase delta interaction protein) (p38), POLDIP2, PDIP38 POLD4
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	POLDIP2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.