



GPSM3 Antibody, FITC conjugated

Product Code	CSB-PA896517LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y4H4
Immunogen	Recombinant Human G-protein-signaling modulator 3 protein (1-160AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	G-protein-signaling modulator 3 (Activator of G-protein signaling 4) (G18.1b) (Protein G18), GPSM3, AGS4 C6orf9 G18
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	GPSM3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.