



TIMM10B Antibody, Biotin conjugated

Product Code	CSB-PA896532LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y5J6
Immunogen	Recombinant Human Mitochondrial import inner membrane translocase subunit Tim10 B protein (1-103AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mitochondrial import inner membrane translocase subunit Tim10 B (Fracture callus protein 1) (FxC1) (Mitochondrial import inner membrane translocase subunit Tim9 B) (TIMM10B) (Tim10b), TIMM10B, FXC1 TIM9B TIMM9B
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	TIMM10B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.