



PPP2R3B Antibody, HRP conjugated

Product Code	CSB-PA896533HB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y5P8
Immunogen	Recombinant Human Serine/threonine-protein phosphatase 2A regulatory subunit B" subunit beta protein (1-176AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serine/threonine-protein phosphatase 2A regulatory subunit B" subunit beta (PP2A subunit B isoform PR48) (Protein phosphatase 2A 48 kDa regulatory subunit), PPP2R3B, PPP2R3L
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PPP2R3B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.