



PPP2R3B Antibody, HRP conjugated

Product Code	CSB-PA896533HB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y5P8
Immunogen	Recombinant Human Serine/threonine-protein phosphatase 2A regulatory subunit B" subunit beta protein (1-176AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serine/threonine-protein phosphatase 2A regulatory subunit B" subunit beta (PP2A subunit B isoform PR48) (Protein phosphatase 2A 48 kDa regulatory subunit), PPP2R3B, PPP2R3L
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PPP2R3B