



TPRKB Antibody, FITC conjugated

Product Code	CSB-PA896710LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y3C4
Immunogen	Recombinant Human EKC/KEOPS complex subunit TPRKB protein (1-175AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	EKC/KEOPS complex subunit TPRKB (PRPK-binding protein) (TP53RK-binding protein), TPRKB
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TPRKB