



COQ4 Antibody, FITC conjugated

Product Code	CSB-PA896885LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y3A0
Immunogen	Recombinant Human Ubiquinone biosynthesis protein COQ4 homolog, mitochondrial protein (31-265AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ubiquinone biosynthesis protein COQ4 homolog, mitochondrial (Coenzyme Q biosynthesis protein 4 homolog), COQ4
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	COQ4