



# MRFAP1 Antibody, Biotin conjugated

Product Code	CSB-PA896920ED01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y605
Immunogen	Recombinant Human MORF4 family-associated protein 1 protein (1-127AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	MORF4 family-associated protein 1 (Protein PGR1) (Protein associated with MRG of 14 kDa), MRFAP1, PAM14 PGR1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	MRFAP1