



MRFAP1 Antibody, Biotin conjugated

Product Code	CSB-PA896920ED01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y605
Immunogen	Recombinant Human MORF4 family-associated protein 1 protein (1-127AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	MORF4 family-associated protein 1 (Protein PGR1) (Protein associated with MRG of 14 kDa), MRFAP1, PAM14 PGR1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	MRFAP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.