




SLC17A4 Antibody

Product Code	CSB-PA897087LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y2C5
Immunogen	Recombinant Human Probable small intestine urate exporter protein (14-57AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IF; Recommended dilution: IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Probable small intestine urate exporter (Solute carrier family 17 member 4), SLC17A4
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC17A4

Image



Immunofluorescent analysis of HeLa cells using CSB-PA897087LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.