



SLC17A4 Antibody, Biotin conjugated

Product Code	CSB-PA897087LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y2C5
Immunogen	Recombinant Human Probable small intestine urate exporter protein (14-57AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Probable small intestine urate exporter (Solute carrier family 17 member 4), SLC17A4
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC17A4