



NDUFAF1 Antibody, FITC conjugated

Product Code	CSB-PA897101LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y375
Immunogen	Recombinant Human Complex I intermediate-associated protein 30, mitochondrial protein (25-327AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Complex I intermediate-associated protein 30, mitochondrial (NADH dehydrogenase [ubiquinone] 1 alpha subcomplex assembly factor 1), NDUFAF1, CIA30
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	NDUFAF1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.