



NDUFAF1 Antibody, Biotin conjugated

Product Code	CSB-PA897101LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y375
Immunogen	Recombinant Human Complex I intermediate-associated protein 30, mitochondrial protein (25-327AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Complex I intermediate-associated protein 30, mitochondrial (NADH dehydrogenase [ubiquinone] 1 alpha subcomplex assembly factor 1), NDUFAF1, CIA30
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	NDUFAF1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.