



NDUFAF1 Antibody, Biotin conjugated

Product Code	CSB-PA897101LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y375
Immunogen	Recombinant Human Complex I intermediate-associated protein 30, mitochondrial protein (25-327AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Complex I intermediate-associated protein 30, mitochondrial (NADH dehydrogenase [ubiquinone] 1 alpha subcomplex assembly factor 1), NDUFAF1, CIA30
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	NDUFAF1