



FARP1 Antibody, Biotin conjugated

Product Code	CSB-PA897492LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y4F1
Immunogen	Recombinant Human FERM, ARHGEF and pleckstrin domain-containing protein 1 protein (1-300AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	FERM, ARHGEF and pleckstrin domain-containing protein 1 (Chondrocyte-derived ezrin-like protein) (FERM, RhoGEF and pleckstrin domain-containing protein 1) (Pleckstrin homology domain-containing family C member 2) (PH domain-containing family C member 2), FARP1, CDEP PLEKHC2
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	FARP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.