




RAB11B Antibody

Product Code	CSB-PA614903LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15907
Immunogen	Recombinant Human Ras-related protein Rab-11B protein (2-215AA)
Raised In	Rabbit
Species Reactivity	Human, Rat
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ras-related protein Rab-11B (GTP-binding protein YPT3), RAB11B, YPT3
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	RAB11B

Image



Western blot
 All lanes: RAB11B antibody at 2µg/ml
 Lane 1: Hela whole cell lysate
 Lane 2: Rat gonadal tissue
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 25, 20 kDa
 Observed band size: 25 kDa



Immunohistochemistry of paraffin-embedded human placenta tissue using CSB-PA614903LA01HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.