



Recombinant Mouse BAG family molecular chaperone regulator 3 (Bag3)

| | |
|--------------------------|---|
| Product Code | CSB-MP864372MO |
| Abbreviation | Bag3 |
| Storage | Store at -20°C, for extended storage, conserve at -20°C or -80°C. |
| Uniprot No. | Q9JLV1 |
| Product Type | Recombinant Protein |
| Immunogen Species | Mus musculus (Mouse) |
| Purity | >85% (SDS-PAGE) |
| Sequence | SAATQSPMMQMASGNGASDRDPLPPGWEIKIDPQTGWPFVVDHNSRTTTWN DPRVPPEGPKDTASSANGPSRDGSROLLPIREGHPIYPQLRPGYIPIVLHEGSE NRQPHLFHAYSQPGVQRFRTAAAATPQRSQSPLRGGMTEAAQTDKQCGQ MPATATTAQAQPPTAHGPERSQSPAASDCSSSSSSASLPSSGRSSLGSHQLP RGYIPIVIHEQNITRPAAQPSFHQAQKTHYPAAQQGEYQPQQPVYHKIQGDDW EPRPLRAASPFSPVRGASSREGSPARSGTPVHCPSPIRVHTVDRPQPMTH REPPVVTQPENKPEKPGPAGPDLPPGHIPIQVIRREADSKPVSQKSPPPAEK VEVKVSSAPIPCSPSPAPSAVPSPPKNVAEEQKAAPSPAPAEPAPKSGEAE TPPKHPGVLKVEAILEKVQGLEQAVDSFEGKKTDDKYLMIIEYLTKEALLDSV DPEGRADVRQARRDGVKRVQVILEKLEKQKIDVPGQVQVYELQPSNLEAEQP LQEIMGAVVADKDKKGPENKDPQTESQQLEAKAATPPNPSNPADSAGNLVAP |
| Research Area | Others |
| Source | Mammalian cell |
| Target Names | Bag3 |
| Protein Names | Recommended name: BAG family molecular chaperone regulator 3 Short name= BAG-3 Alternative name(s): Bcl-2-associated athanogene 3 Bcl-2-binding protein Bis |
| Expression Region | 2-577aa |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full Length of Mature Protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. |



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.