



Recombinant Human Delta-like protein 3 (DLL3)?partial

Product Code	CSB-BP882142HU
Abbreviation	DLL3
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9NYJ7
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	AGVFELQIHSFGPGPGAPRSPCSARLPCRLFFRVCLKPGLSEEAAESPCAL GAALSARGPVYTEQPGAPAPDLPLPDGLLQVPPFRDAWPGTFSFIETWREELG DQIGGPAWSLLARVAGRRRLAAGGPWARDIQRAGAWELRFSYRARCEPPAV GTACTRLCRPRSAPSRGGLRPCAPLEDECEAPLVCRAGCSPEHGFCCEQP GECRCLEGWTGPLCTVPVSTSSCLSPRGPSSATTGCLVPGPGPCDGNPCAN GGSCSETPRSFECTCPRGFYGLRCEVSGVTCADGPCFNGGLCVGGADPDSA YICHCPPGFQGSNCEKRVDRCSLQPCRNGGLCLDLGHALRCRCRAGFAGPR CEHDLDDCAGRACANGGTCVEGGGAHRCSCALGFGGRDRCRERADPCAARP CAHGGRCYAHFSGLVACAPGYMGARCEFPVHPDGASALPAAPPGLRPGDP QRYL
Research Area	Developmental Biology
Source	Baculovirus
Target Names	DLL3
Protein Names	Recommended name: Delta-like protein 3 Alternative name(s): Drosophila Delta homolog 3 Short name= Delta3
Expression Region	27-492aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Extracellular domain
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.