



# Recombinant Human Inactive N-acetylated-alpha-linked acidic dipeptidase-like protein 2 (NAALADL2)

<b>Product Code</b>	CSB-MP704507HU
<b>Abbreviation</b>	NAALADL2
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	Q58DX5
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	<p>           GYYVHTNCPSDAPSSGTVDPQLYQEILKTIQAEDIKKSFRNLVQLYKNEDDMEI            SKKIKTQWTSLGLEDVQFVNYSVLLDLPGPSSTVTLSSSGQCFHPNGQPCSE            EARKDSSQDLLYSYAAYSAGKTLKAEVIDVSYGMADDLKRIRKIKNVNQNIALLK            LGKLPPLYKLSSEKAGFGGVLLYIDPCDLPKTVNPSHDTFMVSLNPGGDPSTP            GYPSVDESEFRQSRSNLTSLLVQPISAPLVAKLISSPKARTKNEACSSLELPNNEI            RVVSMQVQTVTKLKTVTNVVGFVMGLTSPDRYIIVGSHHHTAHSYNGQEWAS            STAITAFIRALMSKVKRGWRPDRTIVFCSWGGTAFGNIGSYEWGEDFKKVLQ            KNVVAYISLHSPIRGNSSLYPVASPSLQQLVVEKNNFNCTRRAQCPETNISSIQI            QGDADYFINHLGVPIVQFAYEDIKTEGSPFLSEARFSTRATKIEEMDPSFNLHE            TITKLSGEVILQIANEPVLPFNALDIALEVQNNLKGDPNTHQLLAMALRLRESA            ELFQSDEMPPANDPKERAPIRIRMLNDILQDMEKSFLVKQAPPGFYRNILYHLD            EKTSRFSILIEAWEHCKPLASNETLQEALSEVLNSINSAQVYFKAGLDVFKSVLD            GKN         </p>
<b>Research Area</b>	Cell Biology
<b>Source</b>	Mammalian cell
<b>Target Names</b>	NAALADL2
<b>Protein Names</b>	Recommended name: Inactive N-acetylated-alpha-linked acidic dipeptidase-like protein 2 Short name= NAALADase L2
<b>Expression Region</b>	143-795aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Extracellular domain
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.