



# Recombinant Influenza B virus Nucleoprotein (NP)

|                          |  |
|--------------------------|--|
| <b>Product Code</b>      | CSB-BP361330IJK  |
| <b>Relevance</b>         | Encapsidates the negative strand viral RNA, protecting it from nucleases. The encapsidated genomic RNA is termed the ribonucleoprotein (RNP) and serves as template for transcription and replication. The RNP needs to be localized in the host nucleus to start an infectious cycle, but is too large to diffuse through the nuclear pore complex. NP comprises at least 2 nuclear localization signals that are responsible for the active RNP import into the nucleus through cellular importin alpha/beta pathway. Later in the infection, nuclear export of RNPs are mediated through viral proteins NEP interacting with M1 which binds nucleoproteins. It is possible that nucleoprotein binds directly host exportin-1/XPO1 and plays an active role in RNPs nuclear export. M1 interaction with RNP seems to hide nucleoprotein's nuclear localization signals. Soon after a virion infects a new cell, M1 dissociates from the RNP under acidification of the virion driven by M2 protein. Dissociation of M1 from RNP unmask nucleoprotein's nuclear localization signals, targeting the RNP to the nucleus. |
| <b>Abbreviation</b>      | Recombinant Influenza B virus NP protein   |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |
| <b>Uniprot No.</b>       | P04665   |
| <b>Product Type</b>      | Recombinant Protein  |
| <b>Immunogen Species</b> | Influenza B virus (strain B/Lee/1940)  |
| <b>Purity</b>            | ≥ 85% as determined by SDS-PAGE.   |
| <b>Sequence</b>          | MSNMDIDSINTGTIDKTPEELTPGTSGATRPIIKPATLAPPSNKRTRNPSPERTT<br>TSSETDIGRKIQKKQTPTEIKKSVYKMOVVVLGFEFYNQMMVKAGLNDDMERNLI<br>QNAQAVERRILLAATDDKKTEYQKKRNARDVKEGKEEIDHNKTGGTFYKMVRD<br>DKTIYFSPIKITFLKEEVKTMKTTMGSDGFSGLNHIMIGHSQMNDVCFQRSKG<br>LKRVLGDPSTLIFAGSTLPRRSGTTGVAIKGGGTLVDEAIRFIGRAMADRGLL<br>RDIKAKTAYEKILLNLKNKCSAPQQKALVDQVIGSRNPGIADIEDLTLLARSMVV<br>VRPSVASKVVLPISIYAKIPQLGFNTEEYSVMGYEAMALYNMATPVSILRMGDD<br>AKDKSQLFFMSCFGAAYEDLRVLSALTGTEFKPRSALKCKGFHVPAAKEQVEG<br>MGAALMSIKLQFWAPMTRSGGNEVSGEGGSGQISCSPVFAVERPIALSKQAV<br>RRMLSMNVEGRDADVKGNLLKMMNDSMAKKTSGNAFIGKMKFQISDKNKVN<br>PIEPIKQTIPNFFFGRDTAEDYDDLDDY   |
| <b>Research Area</b>     | Microbiology   |
| <b>Source</b>            | Baculovirus  |
| <b>Target Names</b>      | NP   |
| <b>Protein Names</b>     | Nucleocapsid protein   |
| <b>Expression Region</b> | 1-560aa  |




**Notes** Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

**Tag Info** N-terminal 10xHis-tagged

**Mol. Weight** 64.3kDa

**Protein Length** Full Length

**Image**



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

**Reconstitution** We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

**Shelf Life** The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.