



Recombinant Rat Ribosyldihydronicotinamide dehydrogenase [quinone] (Nqo2)

Product Code	CSB-EP721470RA-B
Abbreviation	Nqo2
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q6AY80
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	>85% (SDS-PAGE)
Sequence	MAGKKVLLVYAHQEPKSFNGSMKQVAVEELSKQGCTVTVSDLYTMNFEP RAT RNDVTGALSNPEVFKYGIEAYEAYKKKALTS DILEEQRKVQEADLVIFQFPLYWFSVPAILK GWMMDRVLCCQGFADFVPGFYDSGFLKDKL ALLSFTTGGTAEMYTKAGVNGDFRYFLWPL LQHGTLHFCGFKVLAPQISFGPEVSSEEQRK VMLASWVQR LKSIWKEEPIHCTPSWYFQG
Source	E.coli
Target Names	Nqo2
Protein Names	Recommended name: Ribosyldihydronicotinamide dehydrogenase [quinone] EC= 1.10.99.2 Alternative name(s): NRH dehydrogenase [quinone] 2 NRH:quinone oxidoreductase 2 Quinone reductase 2 Short name= QR2
Expression Region	1-231aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.