



Recombinant Human Tyrosine-protein kinase transmembrane receptor ROR1 (ROR1)

Product Code	CSB-MP020067HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q01973
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	QETELSVSAELVPTSSWNISSELNKDSYLTLEPMNNITSLGQTAELHCKVSG NPPPTIRWFKNDAPVVQEPRRLSFRSTIYGSRLRIRNLDTTDTGYFQCVATNG KEVVSSTGVLVFKFGPPPTASPGYSDEYEEDGFCQPYRGIACARFIGNRTVYM ESLHMQGEIENQITAAFTMIGTSSHLSDKCSQFAIPSLCHYAFPYCDETSSVPK PRDLCRDECEILENVLCQTEYIFARSNPMILMRLKLPNCEDLPQPESPEAANCI RIGIPMADPINKNHKCYNSTGVDYRGTVSVTKSGRQCQPWNSQYPHTHTFTA LRFPELNGGHSYCRNPGNQKEAPWCFTLDENFKSDLCDIPAC
Research Area	Neuroscience
Source	Mammalian cell
Target Names	ROR1
Protein Names	Recommended name: Tyrosine-protein kinase transmembrane receptor ROR1 EC= 2.7.10.1Alternative name(s): Neurotrophic tyrosine kinase, receptor-related 1
Expression Region	30-391aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Extracellular domain
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.