



Recombinant Human Collagen alpha-3 (VI) chain (COL6A3), partial

Product Code	CSB-BP005753HU1
Relevance	Collagen VI acts as a cell-binding protein.
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P12111
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	KPMVKMSREVQVFEITENSAKLHWERAEPGPYFYDLTVTSAHDQSLVLKQN LVTDRVIGLLAGQTYHVAVVCYLRSQVRATYHGSFSTKKSQPPPPQPARSA SSSTINLMVSTEPLALTETDICKLPKDEGTCRDFILKWYYDPNTKSCARFWYGG CGGNENKFGSQKECEKVCAPVLAKPGVISVMG
Source	Baculovirus
Target Names	COL6A3
Expression Region	2986-3176aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Partial
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.