




# Recombinant Human Gelsolin (GSN)

<b>Product Code</b>	CSB-EP009965HU
<b>Abbreviation</b>	Recombinant Human GSN protein
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P06396
<b>Form</b>	Liquid or Lyophilized powder
<b>Storage Buffer</b>	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
<b>Product Type</b>	Recombinant Human Gelsolin (GSN)
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Sensitivity</b>	Not Test
<b>Purity</b>	Greater than 95% as determined by SDS-PAGE.
<b>Sequence</b>	<p>           ATASRGASQAGAPQGRVPEARPNMSMVVEHPEFLKAGKEPGLQIWRVEKFDLV            PVPTNLYGDDFFTGDAYVILKTVQLRNGNLQYDLHYWLGNECSQDESGAAAIPT            VQLDDYLNGRAVQHREVVQGFESATFLGYFKSGLKYKKGGSVASGFKHVVPNEV            VVQRLFQVKGRRVVRATEVPVSWESFNNGDCFILDGNNIHWCGSNSNRYE            RLKATQVSKGIRDNERSGRARVHVSEEGTEPEAMLQVLGPKPALPAGTEDTA            KEDANRKLAKLYKVSNGAGTMSVSLVADENPFAQGALKSEDCFILDHGKDG            KIFVWKGKQANTEERKAALKTASDFITKMDYPKQTQVSVLPEGGETPLFKQFF            KNWRDPDQTDGLGLSYLSSHIANVERVPFDAATLHTSTAMAAQHGMDDDG            GQKQIWRIEGSNKVPVDPATYGGFYGGDSYIILYNRHHGGRQGGIYNWQGA            QSTQDEVAASAILTAQLDEELGGTPVQSRVVQGKEPAHLMSLFGGKPMIYK            GTSREGGQTAPASTRLFQVRANSAGATRAVEVLPKAGALNSNDAFVLKTPSA            AYLWVGTGASEAEKTGAQELLRLRAQPVQVAEGSEPDGFWEALGGKAAAYR            TSPRLKDKKMDAHPRLFACSNKIGRFVIEEVPGELMQEDLATDDVMLLDTWD            QVFVWVGKDSQEEEEKTEALTSAKRYIETDPANRDRRTPITVVKQGFEPSPFVG            WFLGWDDDDYWSVDPLDRAMAELAA         </p>
<b>Research Area</b>	Signal Transduction
<b>Source</b>	E.coli
<b>Target Names</b>	GSN
<b>Expression Region</b>	28-782aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 6xHis-tagged
<b>Mol. Weight</b>	87.0 kDa

**Protein Length**

Full Length of Mature Protein

**Image**

(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

**Reconstitution**

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.