






Recombinant Nicotiana tabacum Osmotin (AP24)

Product Code	CSB-EP320333NHE
Abbreviation	Recombinant Nicotiana tabacum AP24 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P14170
Product Type	Recombinant Protein
Immunogen Species	Nicotiana tabacum (Common tobacco)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	ATIEVRNNCPYTVWAASTPIGGGRRRLDRGQTWVINAPRGTKMARVWGRNTC NFNAAGRGTCTGDCGGVLQCTGWGKPPNTLAEYALDQFSGLDLFDISLVD GFNIPMTFAPTNPSSGGKCHAIHCTANINGECPRELVPGGCNPCTTFGGQQ YCCTQGPGPTFFSKFFKQRCPDAYSYPQDDPTSTFTCPGGSTNYRVIFCPN GQAHPNFPLEMPGSDEVAK
Research Area	Others
Source	E.coli
Target Names	AP24
Expression Region	22-246aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-SUMO-tagged
Mol. Weight	40.3kDa
Protein Length	Full Length of Mature Protein

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.