




# Recombinant Human Glucagon-like peptide 1 receptor (GLP1R), partial

|                          |                                                                                                                                                                                                                                                                                 |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>      | CSB-YP009514HU                                                                                                                                                                                                                                                                  |
| <b>Relevance</b>         | This is a receptor for glucagon-like peptide 1. The activity of this receptor is mediated by G proteins which activate adenyl cyclase.                                                                                                                                          |
| <b>Abbreviation</b>      | Recombinant Human GLP1R protein, partial                                                                                                                                                                                                                                        |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| <b>Uniprot No.</b>       | P43220                                                                                                                                                                                                                                                                          |
| <b>Product Type</b>      | Recombinant Protein                                                                                                                                                                                                                                                             |
| <b>Immunogen Species</b> | Homo sapiens (Human)                                                                                                                                                                                                                                                            |
| <b>Purity</b>            | ≥ 90% as determined by SDS-PAGE.                                                                                                                                                                                                                                                |
| <b>Research Area</b>     | Cardiovascular                                                                                                                                                                                                                                                                  |
| <b>Source</b>            | Yeast                                                                                                                                                                                                                                                                           |
| <b>Target Names</b>      | GLP1R                                                                                                                                                                                                                                                                           |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                                                                                                                                                             |
| <b>Tag Info</b>          | N-terminal His-tagged<br>Tag-Free                                                                                                                                                                                                                                               |
| <b>Mol. Weight</b>       | 16.3kDa                                                                                                                                                                                                                                                                         |
| <b>Protein Length</b>    | Partial                                                                                                                                                                                                                                                                         |

## Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

## Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

---

### **Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.