




# Recombinant Human IgG receptor FcRn large subunit p51 (FCGRT), partial

<b>Product Code</b>	CSB-YP008545HU
<b>Relevance</b>	Binds to the Fc region of monomeric immunoglobulins gamma. Mediates the uptake of IgG from milk. Possible role in transfer of immunoglobulin G from mother to fetus.
<b>Abbreviation</b>	Recombinant Human FCGRT protein, partial
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P55899
<b>Alias</b>	IgG Fc fragment receptor transporter alpha chain Neonatal Fc receptor
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥ 90% as determined by SDS-PAGE.
<b>Sequence</b>	AESHLSLLYHLTAVSSPAPGTPAFWVSGWLGPQQYLSYNLSLRGEAEP CGAW VWENQVSWYWEKETTDLRIKEKLFLEAFKALGGKGPYTLQGLLGCELGPDNT SVPTAKFALNGEEFMNFDLKQGTWGGDWPEALAISQRWQQQDKAANKELTF LLFSCPHRLREHLERGRGNLEWKEPPSMRLKARPSSPGFSVLTCSAFSFYPP ELQLRFLRNGLAAGTGQGDGFGPNSDGSFHASSSLTVKSGDEHHYCCIVQHAG LAQPLRVELESPAKSS
<b>Research Area</b>	Immunology
<b>Source</b>	Yeast
<b>Target Names</b>	FCGRT
<b>Protein Names</b>	Recommended name: IgG receptor FcRn large subunit p51 Short name= FcRn Alternative name(s): IgG Fc fragment receptor transporter alpha chain Neonatal Fc receptor
<b>Expression Region</b>	24-297aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 6xHis-tagged
<b>Mol. Weight</b>	32.4kDa
<b>Protein Length</b>	Extracellular Domain
<b>Image</b>	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

### Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.