



L1 Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA08349A0Rb |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P03101 |
| Immunogen | Recombinant Human papillomavirus type 16 Major capsid protein L1 protein (6-471AA) |
| Raised In | Rabbit |
| Species Reactivity | Human papillomavirus type 16 |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Major capsid protein L1, L1 |
| Immunogen Species | Human papillomavirus type 16 |
| Research Area | Microbiology |
| Target Names | L1 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |