



PPIA Antibody Pair

Product Code	CSB-EAP07814
Species Reactivity	Human,Rat,Mouse,Guinea pig,Bovine,Pig,Dog,Sheep
Tested Applications	S-ELISA
Form	Liquid
Alias	CYPA,PPIA,Cyclophilin A Cyclosporin A-binding protein Rotamase A
Product Type	Antibody Pairs
Immunogen Species	Homo sapiens (Human)
Gene Names	PPIA
Protein Names	Peptidyl-prolyl cis-trans isomerase A
Notes	We recommend using the capture antibody at a concentration of 1.5ug/ml and the detection antibody at a concentration of 0.67ug/ml.Optimal dilutions should be determined experimentally by the researcher.
Host	Capture: Rabbit Detection: Rabbit
Components	Capture: CSB-EAP07814C Detection: CSB-EAP07814D(Biotin) Reagents are sufficient for at least 5 x 96 well plates using recommended protocol.
Storage-Buffer	Capture: 50% Glycerol, 0.01M PBS, PH 7.4 Detection: 50% Glycerol, 0.01M PBS, PH 7.4