



YWHAB Antibody Pair

Product Code	CSB-EAP01544
Species Reactivity	Human,Rat,Mouse,Guinea pig,Bovine,Horse,Pig,Dog,Chicken,Goat,Sheep
Tested Applications	S-ELISA
Form	Liquid
Alias	Protein 1054,Protein kinase C inhibitor protein 1,KCIP-1,14-3-3 protein beta/alpha, N-terminally processed
Product Type	Antibody Pairs
Immunogen Species	Homo sapiens (Human)
Target Names	YWHAB
Protein Names	14-3-3 protein beta/alpha
Notes	We recommend using the capture antibody at a concentration of 5ug/ml and the detection antibody at a concentration of 1.25ug/ml.Optimal dilutions should be determined experimentally by the researcher.
Host	Capture: Rabbit Detection: Rabbit
Components	Capture: CSB-EAP01544C Detection: CSB-EAP01544D(Biotin) Reagents are sufficient for at least 5 x 96 well plates using recommended protocol.
Storage-Buffer	Capture: 50% Glycerol, 0.01M PBS, PH 7.4 Detection: 50% Glycerol, 0.01M PBS, PH 7.4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.