



CTNNB1 Antibody Pair

Product Code	CSB-EAP00174
Species Reactivity	Human,Rat,Mouse,Guinea pig,Bovine,Goat,Sheep
Tested Applications	S-ELISA
Form	Liquid
Alias	CTNNB1, CTNNB, OK/SW-cl.35, PRO2286,Beta-catenin
Product Type	Antibody Pairs
Immunogen Species	Homo sapiens (Human)
Gene Names	CTNNB1
Protein Names	Catenin beta-1
Notes	We recommend using the capture antibody at a concentration of 1ug/ml and the detection antibody at a concentration of 0.1ug/ml.Optimal dilutions should be determined experimentally by the researcher.
Host	Capture: Rabbit Detection: Rabbit
Components	Capture: CSB-EAP00174C Detection: CSB-EAP00174D(Biotin) Reagents are sufficient for at least 5 x 96 well plates using recommended protocol.
Storage-Buffer	Capture: 50% Glycerol, 0.01M PBS, PH 7.4 Detection: 50% Glycerol, 0.01M PBS, PH 7.4