



CSTB Antibody Pair

Product Code	CSB-EAP00224
Species Reactivity	Human,Mouse,Pig,Dog,Rabbit,Goat,Duck
Tested Applications	S-ELISA
Form	Liquid
Alias	CST6, STFB,CSTB,CPI-B,Liver thiol proteinase inhibitor Stefin-B
Product Type	Antibody Pairs
Immunogen Species	Homo sapiens (Human)
Gene Names	CSTB
Protein Names	Cystatin-B
Notes	We recommend using the capture antibody at a concentration of 0.2ug/ml and the detection antibody at a concentration of 0.1ug/ml.Optimal dilutions should be determined experimentally by the researcher.
Host	Capture: Rabbit Detection: Rabbit
Components	Capture: CSB-EAP00224C Detection: CSB-EAP00224D(Biotin) Reagents are sufficient for at least 5 x 96 well plates using recommended protocol.
Storage-Buffer	Capture: 50% Glycerol, 0.01M PBS, PH 7.4 Detection: 50% Glycerol, 0.01M PBS, PH 7.4