



CP Antibody Pair

Product Code	CSB-EAP14599
Species Reactivity	Human,Hamster,Mouse,Bovine,Horse,Pig,Dog,Rabbit
Tested Applications	S-ELISA
Form	Liquid
Alias	Ferroxidase, CP
Product Type	Antibody Pairs
Immunogen Species	Homo sapiens (Human)
Gene Names	CP
Protein Names	Ceruloplasmin
Notes	We recommend using the capture antibody at a concentration of 2ug/ml and the detection antibody at a concentration of 0.25ug/ml.Optimal dilutions should be determined experimentally by the researcher.
Host	Capture: Rabbit Detection: Rabbit
Components	Capture: CSB-EAP14599C Detection: CSB-EAP14599D(Biotin) Reagents are sufficient for at least 5 x 96 well plates using recommended protocol.
Storage-Buffer	Capture: 50% Glycerol, 0.01M PBS, PH 7.4 Detection: 50% Glycerol, 0.01M PBS, PH 7.4