



DAP Antibody Pair

Product Code	CSB-EAP01034
Species Reactivity	Human,Hamster,Pig,Dog,Rabbit
Tested Applications	S-ELISA
Form	Liquid
Alias	DAP-1
Product Type	Antibody Pairs
Immunogen Species	Homo sapiens (Human)
Gene Names	DAP
Protein Names	Death-associated protein 1
Notes	We recommend using the capture antibody at a concentration of 0.8ug/ml and the detection antibody at a concentration of 0.4ug/ml.Optimal dilutions should be determined experimentally by the researcher.
Host	Capture: Rabbit Detection: Rabbit
Components	Capture: CSB-EAP01034C Detection: CSB-EAP01034D(Biotin) Reagents are sufficient for at least 5 x 96 well plates using recommended protocol.
Storage-Buffer	Capture: 50% Glycerol, 0.01M PBS, PH 7.4 Detection: 50% Glycerol, 0.01M PBS, PH 7.4