



IL8 Antibody Pair

Product Code	CSB-EAP08327
Species Reactivity	Human,Rat,Mouse,Guinea pig,Bovine,Pig,Dog,Chicken,Goat,Sheep
Tested Applications	S-ELISA
Form	Liquid
Alias	CXCL8,IL8 IL-8 C-X-C motif chemokine 8 Emotakin Granulocyte chemotactic protein 1 GCP-1 Monocyte-derived neutrophil chemotactic factor MDNCF Monocyte-derived neutroph
Product Type	Antibody Pairs
Immunogen Species	Homo sapiens (Human)
Gene Names	CXCL8
Protein Names	Interleukin-8
Notes	We recommend using the capture antibody at a concentration of 0.5ug/ml and the detection antibody at a concentration of 0.25ug/ml.Optimal dilutions should be determined experimentally by the researcher.
Host	Capture: Rabbit Detection: Rabbit
Components	Capture: CSB-EAP08327C Detection: CSB-EAP08327D(Biotin) Reagents are sufficient for at least 5 x 96 well plates using recommended protocol.
Storage-Buffer	Capture: 50% Glycerol, 0.01M PBS, PH 7.4 Detection: 50% Glycerol, 0.01M PBS, PH 7.4