



ORM1 Antibody Pair

Product Code	CSB-EAP05709
Species Reactivity	Human,Rat,Mouse,Guinea pig,Bovine,Horse,Pig,Dog,Rabbit,Sheep
Tested Applications	S-ELISA
Form	Liquid
Alias	AGP1,ORM1,Orosomucoid-1,OMD 1,AGP 1
Product Type	Antibody Pairs
Immunogen Species	Homo sapiens (Human)
Gene Names	ORM1
Protein Names	Alpha-1-acid glycoprotein 1
Notes	We recommend using the capture antibody at a concentration of 0.6ug/ml and the detection antibody at a concentration of 0.34ug/ml.Optimal dilutions should be determined experimentally by the researcher.
Host	Capture: Rabbit Detection: Rabbit
Components	Capture: CSB-EAP05709C Detection: CSB-EAP05709D(Biotin) Reagents are sufficient for at least 5 x 96 well plates using recommended protocol.
Storage-Buffer	Capture: 50% Glycerol, 0.01M PBS, PH 7.4 Detection: 50% Glycerol, 0.01M PBS, PH 7.4