



IL1B Antibody Pair

Product Code	CSB-EAP06574
Immunogen	Rabbit
Species Reactivity	Human?Rat?Guinea pig?Pig?Dog?Rabbit?Goat
Tested Applications	S-ELISA
Form	Liquid
Product Type	Antibody Pairs
Immunogen Species	Homo sapiens (Human)
Protein Names	Interleukin-1 beta
Notes	We recommend using the capture antibody at a concentration of 0.5ug/ml and the detection antibody at a concentration of 1ug/ml.Optimal dilutions should be determined experimentally by the researcher.
Host	Capture: Rabbit Detection: Rabbit
Components	Capture: CSB-EAP06574C Detection: CSB-EAP06574D(Biotin) Reagents are sufficient for at least 5 x 96 well plates using recommended protocol.
Storage-Buffer	Capture: 50% Glycerol, 0.01M PBS, PH 7.4 Detection: 50% Glycerol, 0.01M PBS, PH 7.4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.