



C21orf33 Antibody Pair

Product Code	CSB-EAP01814
Species Reactivity	Human,Mouse,Pig,Dog,Rabbit
Tested Applications	S-ELISA
Form	Liquid
Alias	Protein GT335 Protein KNP-I
Product Type	Antibody Pairs
Immunogen Species	Homo sapiens (Human)
Gene Names	C21orf33
Protein Names	ES1 protein homolog, mitochondrial
Notes	We recommend using the capture antibody at a concentration of 0.5ug/ml and the detection antibody at a concentration of 0.34ug/ml.Optimal dilutions should be determined experimentally by the researcher.
Host	Capture: Rabbit Detection: Rabbit
Components	Capture: CSB-EAP01814C Detection: CSB-EAP01814D(Biotin) Reagents are sufficient for at least 5 x 96 well plates using recommended protocol.
Storage-Buffer	Capture: 50% Glycerol, 0.01M PBS, PH 7.4 Detection: 50% Glycerol, 0.01M PBS, PH 7.4