



RPL17 Antibody Pair

Product Code	CSB-EAP02504
Species Reactivity	Human, Guinea pig, Pig, Dog, Rabbit
Tested Applications	S-ELISA
Form	Liquid
Alias	60S ribosomal protein L23, PD-1, RPL17
Product Type	Antibody Pairs
Immunogen Species	Homo sapiens (Human)
Target Names	RPL17
Protein Names	60S ribosomal protein L17
Notes	We recommend using the capture antibody at a concentration of 0.5ug/ml and the detection antibody at a concentration of 0.125ug/ml. Optimal dilutions should be determined experimentally by the researcher.
Host	Capture: Rabbit Detection: Rabbit
Components	Capture: CSB-EAP02504C Detection: CSB-EAP02504D(Biotin) Reagents are sufficient for at least 5 x 96 well plates using recommended protocol.
Storage-Buffer	Capture: 50% Glycerol, 0.01M PBS, PH 7.4 Detection: 50% Glycerol, 0.01M PBS, PH 7.4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.