



ERP29 Antibody Pair

Product Code	CSB-EAP05345
Species Reactivity	Human,Guinea pig,Bovine,Pig,Dog,Sheep
Tested Applications	S-ELISA
Form	Liquid
Alias	ERp29,Endoplasmic reticulum resident protein 28,ERp28,Endoplasmic reticulum resident protein 31,ERp31,ERP29,C12orf8, ERP28
Product Type	Antibody Pairs
Immunogen Species	Homo sapiens (Human)
Target Names	ERP29
Protein Names	Endoplasmic reticulum resident protein 29
Notes	We recommend using the capture antibody at a concentration of 1ug/ml and the detection antibody at a concentration of 0.5ug/ml.Optimal dilutions should be determined experimentally by the researcher.
Host	Capture: Rabbit Detection: Rabbit
Components	Capture: CSB-EAP05345C Detection: CSB-EAP05345D(Biotin) Reagents are sufficient for at least 5 x 96 well plates using recommended protocol.
Storage-Buffer	Capture: 50% Glycerol, 0.01M PBS, PH 7.4 Detection: 50% Glycerol, 0.01M PBS, PH 7.4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.