



H3F3A Antibody Pair

Product Code	CSB-EAP01435
Species Reactivity	Human,Rat,Mouse,Guinea pig,Bovine,Pig,Dog,Rabbit,Goat,Sheep
Tested Applications	S-ELISA
Form	Liquid
Alias	H3F3B,H3.3B,H3F3A,H3.3A, H3F3,PP781,Histone H3.3
Product Type	Antibody Pairs
Immunogen Species	Homo sapiens (Human)
Target Names	H3F3A
Protein Names	Histone H3.3
Notes	We recommend using the capture antibody at a concentration of 0.3ug/ml and the detection antibody at a concentration of 0.4ug/ml.Optimal dilutions should be determined experimentally by the researcher.
Host	Capture: Rabbit Detection: Rabbit
Components	Capture: CSB-EAP01435C Detection: CSB-EAP01435D(Biotin) Reagents are sufficient for at least 5 x 96 well plates using recommended protocol.
Storage-Buffer	Capture: 50% Glycerol, 0.01M PBS, PH 7.4 Detection: 50% Glycerol, 0.01M PBS, PH 7.4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.